hCG Easy (25 mIU/mL)

Intended Use

Alere™ hCG Easy is intended for the qualitative detection of human chorionic gonadotropin (hCG) in urine as an aid in the diagnosis of pregnancy. For professional *in vitro* diagnostic use only.

Introduction

hCG is a glycoprotein hormone produced by the blastocyst.^{1,2} hCG normally begins to be detected in the urine from 7 days after conception. The sudden rapid rise in concentration of hCG in urine following conception makes it an excellent marker for pregnancy.^{3,4}

Test Principle

Urine is added to the absorbent sampler, which contains blue beads attached to anti-hCG monoclonal antibodies. Urine mobilises the labelled antibodies and moves up the test strip, which contains regions of immobilised antibody. If hCG is present at ≥ 25 mIU/mL, a blue Result Line (R) should become visible in the Result Window (see Figure 1). No Result Line (R) indicates the test is negative. A blue Control Line (C) appears in the Result Window to show that the test has worked correctly.

Kit Contents and Storage

20 individually foil wrapped devices. Store at 2-30 °C. Do not use after the expiry date.

Precautions

1. Do not open the foil pouch until ready to test.
2. Do not use devices that have become wet, or if the foil pouch has been damaged.
3. Properly dispose of all contaminated waste such as used tests.
4. Do not use beyond expiration date printed on the outside of the kit carton.

Sample Collection and Storage

A urine sample collected any time of day is suitable but a first morning urine sample is recommended.⁵ Urine samples must be collected in **clean**, dry plastic or glass containers. Samples may be stored in the refrigerator (2-8 °C) for up to 72 hours, or frozen once at -20 °C for up to 3 months. Samples **must** be allowed to reach room temperature (18-30 °C) prior to testing. No sample preparation is required, although particulate matter should be allowed to settle before testing. **Alere™ hCG Easy** is not validated for use with samples containing preservatives.


Assay Procedure

Ensure all **Alere™ hCG Easy** devices and samples are at 18-30 °C. When ready to test, tear open the foil wrapper and remove the device.

After use of device in test, place cap over absorbent sampler.


Do not invert the device during testing

Follow **one** of the following procedures:


Method A – Dip

- Dip the absorbent sampler into the sample to point X, as shown in the diagram.
- **DO NOT** immerse any plastic parts in urine.
- Hold in place for 15 seconds.
- Remove the device from the sample.
- Place the cap over the absorbent sampler.
- **DO NOT** invert the device.
- Read result after 3 minutes. **DO NOT** interpret after 10 minutes.

Method B – Dip And Leave

- Dip **half** of the absorbent sampler into the sample to point Y, as shown in the diagram.
- Leave in place for 3 to 10 minutes.
- Remove the device from the sample.
- Place the cap over the absorbent sampler.
- **DO NOT** invert the device.
- Read result straight away. **DO NOT** interpret after 10 minutes.


Method C – Pipette

- Place the device on a clean, dry, flat surface with the absorbent side of the sampler facing upwards.
- Pipette 100 μ L of the urine sample onto the absorbent sampler.
- Place the cap over the absorbent sampler.
DO NOT invert the device.
- Read result after 3 minutes. **DO NOT** interpret after 10 minutes.

Method D – Urine Stream

- Ask the patient to hold the absorbent sampler pointing downwards in their urine stream for 3 to 7 seconds **only**.
- Place the cap over the absorbent sampler.
DO NOT invert the device.
- Place the device on a flat surface.
- Read result after 3 minutes. **DO NOT** interpret after 10 minutes.

Interpretation of Results

The test can be read 3-10 minutes after applying the sample, regardless of which test procedure is used. See Figure 1.

- A blue Control Line (C) must appear within 3 minutes to show that the test has worked correctly.
- A **positive result**: two blue lines – a Result Line (R) and a Control Line (C) – in the Result Window within 3 minutes. One line may be fainter or darker than the other but this does not affect the interpretation of the result.
- A **negative result**: one blue line – a Control Line (C) – in the Result Window within 3 minutes.
- A **negative result**: one blue line – a Control Line (C) – in the Result Window within 3 minutes.
- An **invalid result**: no blue line in the Control Line (C) region in the Result Window within 3 minutes.
- If no Control Line appears within 3 minutes, the test is invalid and must be repeated with a new **Alere™ hCG Easy** device.
- Any result that appears after 10 minutes must be ignored.


Figure 1

Limitations of the Test

1. Positive results from very early pregnancy may later prove negative as natural termination is estimated to occur in 31% of all conceptions.⁶ It is therefore recommended that weak positive results are re-tested 48-72 hours later with a first morning urine sample.
2. A negative result may be obtained if a urine sample is too dilute. If pregnancy is still suspected, it is recommended the patient should be retested 48-72 hours later with a first morning urine sample.
3. Concentrations of hCG are generally lower in ectopic pregnancy than expected normal values for a given gestational age. Abnormal pregnancy cannot be distinguished from normal pregnancy by hCG levels alone.^{2,7}
4. hCG remains elevated for a time after pregnancy.⁸ Pregnancy tests carried out less than 3 weeks after giving birth or 9 weeks after natural loss or termination may need further evaluation.
5. A number of conditions other than pregnancy can cause elevated levels of urinary hCG e.g. menopause, trophoblastic disease, and certain non-trophoblastic neoplasms.⁹
6. Occasionally, samples containing <25 mIU/mL hCG may test positive but samples containing <5 mIU/mL should be negative. **Alere™ hCG Easy** has been shown to be over 99% accurate.
7. Drugs containing hCG may interfere with **Alere™ hCG Easy**, and produce misleading results.
8. False positive and false negative pregnancy tests may be observed in patients with abnormal bladder or kidney function e.g. enterocystoplasties¹⁰ and renal failure.
9. If the test result with **Alere™ hCG Easy** is not consistent with clinical evidence, further evaluation may be required.
10. Inconsistent results may be obtained if the urine sample contains excessive amounts of bacteria.

Expected Values

Urine samples from healthy males and post-menopausal females generally contain <10 mIU/mL hCG.¹¹ On the first day of the first missed period, the levels of maternal urinary hCG are normally 50-250 mIU/mL. During the first trimester hCG levels peak at up to 200,000 mIU/mL in a typical pregnancy.^{3,12}

Performance Characteristics

Sensitivity

Alere™ hCG Easy can detect hCG in urine at concentrations of ≥25 mIU/mL. This sensitivity has been determined against the 4th International hCG Standard (WHO). Samples containing <5 mIU/mL hCG should be negative.

Prozone Effect

Alere™ hCG Easy has been shown to produce positive results with samples containing up to and including 250,000 mIU/mL hCG, which is higher than the maximum level expected during a typical pregnancy.

Specificity

Alere™ hCG Easy has been evaluated for its cross-reactivity with a variety of substances, including other hormones found in urine. No cross-reactivity was detected when the following substances were added to both positive and negative urine samples:

- | | | |
|------------------------------|------------------------------|---|
| • LH
(500 mIU/mL) | • Bilirubin
(10 µg/mL) | • Salicylic acid
(0.2 mg/mL) |
| • FSH
(200 mIU/mL) | • Ascorbic acid
(2 mg/mL) | • Pregnanediol-3-
glucuronide
(1 µg/mL) |
| • TSH
(1000 µIU/mL) | • Glucose
(10 mg/mL) | • Oestrone-3-
glucuronide
(1 µg/mL) |
| • Haemoglobin
(3.6 mg/mL) | • Uric acid
(0.1 mg/mL) | |
| • Albumin
(12 mg/mL) | • Caffeine
(0.2 mg/mL) | |

Accuracy

In a study conducted using 99 urine samples collected from women for the purposes of pregnancy testing, **Alere™ hCG Easy** was shown to be over 99% accurate and relative specificity and sensitivity were found to be > 99.9%. The study results are summarised below:

[hCG] mIU/mL*	Alere™ hCG Easy	
	+	-
< 5	0	71
5 – 24.9	5	0
≥ 25	23	0

*Assessed by quantitative enzyme immunoassay.

Advice Line

For further information, please contact your distributor, or call Alere Technical Specialists:

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E-post: LAPproductsupport@alere.com

Kanada: Tel: +1 800 818 8335

E-post: CANproductsupport@alere.com

**References / Referencer / Literaturhinweise / Βιβλιογραφία / Referencias bibliográficas /
Lähdeluettelo / Bibliographie / Bibliografia / Literatuur / Referanser / Referências / Referenser**

1. Hsu M., Kolm, P., Leete, J., Dong, K., Muasher, S., Oehninger, S. (1998). Analysis of Implantation in Assisted Reproduction through the use of Serial Human Chorionic Gonadotropin Measurements. *J. Assist. Reprod. Genet.* 15 (8): 496-503.
2. Catt K.J. et al (1975). 40:537-540.
3. Chard T. (1992). Pregnancy tests: review. *Hum. Reprod.* 7 (5): 701-710.
4. Corson S.L. et al (1981). *J. Reprod. Med.*; 26: 611-614.
5. Kaplan L.A. and Pesce A.J. (1989). The C. V. Mosby Company.
6. Wilcox A.J. et al (1988). *N. Engl. J. Med.* 319: 189-194.
7. Braunstein G.D. et al (1978). *Am. J. Obst. Gynec.* 131: 25-32.
8. Steier J. & Myking, O. (1984). Human Chorionic Gonadotropin in Maternal Plasma after Induced Abortion, Spontaneous Abortion, and Removed Ectopic Pregnancy. *Obstet. Gynecol.* 64 (3): 391-394.
9. Braunstein G., Vaitukaitis, J., Carbone, P., Ross, G. (1973). Ectopic Production of Human Chorionic Gonadotrophin by Neoplasms. *Ann. Intern. Med.* 78: 39-45.
10. Nethercliffe J., Trewick, A., Samuel, C., Leaver, R., Woodhouse, C. (2001). False-positive pregnancy tests in patients with enterocystoplasties. *BJU International.* 87 (9) : 780-782.
11. Alfthan H., Haglund, C., Dabek, J., Stenman, U. (1992). Concentrations of Human Choriogonadotropin, its beta-subunit and the core fragment of the beta-subunit in serum and urine of men and nonpregnant women. *Clin. Chem.* 38 (10) : 1981-1987.
12. Lenton, E., Neal, L., Sulaiman, R. (1982). Plasma concentrations of human gonadotropin from the time of implantation until the second week of pregnancy. *Fertil. Steril.* 37 (6) : 773-778.

REF UPK009A

Alere

532 Seventeen Mile Rocks Road
Sinnamon Park, QLD 4073
Australia

www.alere.com

Tel: +61 (0) 7 3363 7100

Fax: +61 (0) 7 3363 7199

MediMark Europe

11, rue Emile Zola – BP 2332
38033 Grenoble Cedex 2 – France

Tel: +33 (0) 476 864 322

Fax: +33 (0) 476 171 982

© 2012 Alere. All rights reserved.

The Alere Logo and Alere are trademarks of the Alere group of companies.

502574/D 2012/06

Attention: By approving the enclosed design draft, you (ABON's Customer) accept all responsibility for the accuracy of the design. If an error is detected following the printing or manufacturing of a material, you (ABON's Customer) are responsible for the cost of any inventory which is deemed unsuitable for sale.

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