FOOD

HOW TO PREPARE

GMP FLAVOUR STANDARDS

Knowing where to begin with flavour training can be a challenge, but we're here to make it simple. Our How to Guide shows you exactly the measurements to use and where to start, helping you build confidence and consistency in every tasting session, whether you're just beginning your sensory journey or are more advanced in flavour evaluation.




BEGINNER	INTERMEDIATE	ADVANCED
1x Capsule per 1L of Water	1x Capsule per 50G of Product	1x Capsule per 100G of Product

FlavorActiV™

A FlavorActiV™ USER GUIDE

How to Prepare GMP Flavour Standards in Food



1. Measure the recommended quantity of final product into a beaker/container.
2. Remove the single capsule needed for the session from its blister pack
3. Hold the capsule upright and at arm's length. Tap the capsule and pull apart each half.* Ensure the capsule is facing away from your face when you pull the capsule apart.
4. Pour the powder into the beaker/container of product.
5. Using a non-odours spoon, mix the flavour standard into the product.
6. Dispense the prepared sample into tasting cups and cover with lid.

**If your capsule is difficult to pull apart- carefully use scissors to cut open one end of the capsule and then pour the solution into the beverage.*

FlavorActiV™